Medtronic

Monopolar Devices

70-10-1352 Rev E
Eff. Date: May 2012

REF 12-201-1 **Endo FB3.0™ Floating Ball**
REF 13-201-1 **Endo SH2.0™ Sealing Hook**

REF 13-121-1 **DS3.5-C™ Dissecting Sealer**
REF 13-101-1 **DS3.0™ Dissecting Sealer**

Follow Instructions for Use

Do Not Reuse

Do Not Restерilize

Humidity Limitation
Operate at 15% – 85%
Non-condensing Humidity

Sterilized by Ethylene Oxide

Does not contain
natural rubber latex

Contains Phthalates

Use by Date

Temperature Limitation
Operate at 50 °F to 86 °F
(10 °C to 30 °C)

Non-Pyrogenic

Medtronic

Manufacturer:

Medtronic Advanced Energy LLC

180 International Drive
Portsmouth, NH 03801 USA

Tel 866.777.9400
Fax 866.222.0900

www.medtronic.com/advancedenergy

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands

Medtronic Australasia Pty Ltd.
97 Waterloo Road
North Ryde NSW 2113 Australia


© Copyright 2011 Medtronic Advanced Energy LLC. All rights reserved. Printed in USA. DS3.0, DS3.5-C, Endo SH2.0 and Endo FB3.0 are registered trademarks of Medtronic Advanced Energy LLC.

ENGLISH


Device Description

A sterile, single-use monopolar device that employs radio-frequency (RF) energy and saline for hemostatic sealing, coagulation and, with the exception of the Endo FB3.0, blunt dissection of soft tissue. The device is equipped with an electrode tip. There is an on/off RF energy actuation button on the handpiece, and saline and electrical lines exit the end of the handpiece.

Figure 1. Side view of an example of a monopolar device.

System Diagram Below is a simple schematic of how the device connects to other instruments. Use the device with any saline delivery system that will deliver 4–20 ml/minute (2–6 ml/minute for Endo FB3.0) of sterile saline (0.9% NaCl) flow to the device. The saline delivery system can be any infusion pump with flow rate setting capability or an intravenous saline bag and intravenous tubing with a 15 drops/ml drip chamber, which can be connected to the device.

Note: Neutral Electrode (NE) = Patient Return Electrode (PRE)

Endo FB3.0 Floating Ball (12-201-1):

Indications For Use

The Endo FB3.0™ device is a sterile, single-use electrosurgical device intended to be used in conjunction with an electrosurgical generator for delivery of radio-frequency energy and saline for hemostatic sealing and coagulation of soft tissue at the operative site. It is intended for, but not limited to, abdominal and thoracic surgery, laparoscopic procedures, endoscopic procedures, and thoracoscopic procedures. The device is not intended for contraceptive tubal coagulation (permanent female sterilization).

Endo SH2.0 Sealing Hook (13-201-1):

Indications For Use

The Endo SH2.0™ sealing hook is a sterile, single-use monopolar electrosurgical device intended to be used in conjunction with an electrosurgical generator for delivery of radio-frequency energy and saline for blunt dissection, hemostatic sealing and coagulation of soft tissue at the operative site. It is intended for, but not limited to, abdominal and thoracic surgery, laparoscopic procedures, and thoracoscopic procedures. The device is not intended for contraceptive tubal coagulation (permanent female sterilization).

DS3.0 and DS3.5-C Dissecting Sealers (13-101-1 and 13-121-1):

Indications For Use

The dissecting sealer is a sterile, single-use monopolar electrosurgical device intended to be used in conjunction with an electrosurgical generator for delivery of radio-frequency energy and saline for blunt dissection, hemostatic sealing and coagulation of soft tissue at the operative site. It is intended for, but not limited to, abdominal and thoracic surgery. The device is not intended for contraceptive tubal coagulation (permanent female sterilization).

Warnings for these devices are:

- The device is not intended for neurosurgical applications.
- Do not activate the device unless it is in contact with tissue to be treated.
- Neuromuscular stimulation can occur with modes which produce electrical arcs between the active electrode and tissue.
- Unless compatible monitoring NE are used with a contact quality monitor, loss of safe contact between the NE and the patient will not result in an auditory alarm.
- Special care should be taken when using the device in the proximity of neural tissue to avoid damage to nerves & similar sensitive structures.
- DO NOT TREAT the following areas at any time:
 - Dura
 - Nerve roots
 - Vertebral end plates after discectomy
 - Skin and skin edges
 - Intact nerves
 - Intact tendons and ligaments
 - Epidural veins
 - Bone surfaces that are intended to be fused
 - Brain tissue
 - Subcutaneous tissue
 - Bone to be covered by implant
- The device is not intended for contraceptive tubal coagulation (permanent female sterilization).
- Ensure the electrode is in contact with the tissue to be treated. Activation and saline flow occur simultaneously. Ensure that saline is flowing at the time of activation. Activating the device with the tip pointing upward may result in inadequate saline flow to the surface intended for treatment.
- **If saline flow stops during the electrosurgical procedure, stop using the device and attempt to resume saline flow. Ensure that the saline source is adequate and the saline delivery system is functioning properly. If unable to resume saline flow, discontinue use and obtain a new device and return the used device to Medtronic Advanced Energy LLC.**
- This device is not intended to be bent, or used as a pry, or any other use not cleared by FDA. Bending or using the device could cause part breakage.
- The device is provided as a sterile, non-pyrogenic, single-use disposable device. **Do not resterilize or reuse this device. Reprocessing, (resterilizing or reusing) this device can result in occluded saline apertures, reducing or preventing flow of saline.**
- Be aware that all exposed metal on the electrode is capable of treating tissue. Use caution to avoid inadvertent treatment of tissue and adjacent structures.
- Be aware that monopolar devices employ RF energy coupled with saline. This coupling effect may result in a deeper tissue effect than conventional RF energy and has the potential for hot saline run off onto delicate structures.
- Protect delicate structures from the hot saline run off by utilization of suction or other protective measures.
- Use suction to minimize the potential for activation of the device in a pool of blood. Activating in a pool of blood may limit the hemostatic effectiveness of the device or increase the risk of an electrode becoming clogged by coagulated blood.
- The tip of the suction wand should not touch the electrode, since this might interfere with the proper function of the device. However, the tip of the suction wand can be as close as 1-2 mm from the electrode when the device is activated.
- Overuse or excessive application of this device may result in contraction, inflammation or necrosis of tissue.
- Use of this device on skin may result in incisional complications such as necrosis or desiccation of the skin.
- Inspect the device and cord for breaks, cracks, nicks, or other damage before use. After connection ensure that they are functioning as intended. Failure to observe this caution may result in injury or electrical shock to the patient or surgical team.
- Interference produced by the operation of the RF surgical equipment may adversely influence the operation of other electronic equipment.
- Failure of the RF surgical equipment could result in an unintended increase of output power.
- Use the lowest setting possible to achieve the desired tissue effect to avoid overtreatment, which could result in swelling, fluid, seroma or unintended tissue necrosis.

Precautions

- In the event that a higher than normal electrosurgical power setting is required, before changing power settings, check all instrument connections, cables and patient contacts. If all connections, cables and patient contacts are fault-free, then increase power settings in small increments, checking carefully after each change.
- Use of a return electrode pad is required for use of this device. Assume that all connections are secure and monitor closely to maintain appropriate patient contact during use. Consult return electrode pad manufacturer's Instructions For Use for proper placement and use.
- **Caution:** When used in solid organ resections, Medtronic monopolar devices are often applied with higher current and relatively long activation times, different from standard, traditional electrosurgical usage. Prior to using higher current, long-duration systems such as Medtronic, assure that all components of the system set-up are compatible, including Medtronic devices, generator, and patient return electrodes (NE). Assure that the area of the NE placement has adequate surface area, musculature, and vasculature for the anticipated current and duration of use. It is important that there be sufficient musculature and vasculature for current distribution so that there are no current concentrations, such as at the edges. Refer to manufacturer's Instructions For Use for all devices used, including the generator, patient return electrode, and all other associated equipment. As stated in the manufacturer's Instructions For Use for NE "Do not place the patient return electrode on anatomy insufficient to allow proper return."

For example, for infants, the safest NE position is on larger anatomy such as the infant's back. Placing a NE around the thigh, calf, or arm increases the possibility of increased temperature at the NE site. **It is also important not to have warming pads or blankets in operation over or around the NE during the time that a Medtronic device is activated** as this could increase the temperature at the site of the NE.

When multiple NE must be used, each NE must be placed in an area that has adequate surface area, musculature, and vasculature. This is important, particularly on patients with small body weight or mass, such as neonatal, infant and pediatric patients, or debilitated adult patients. See NE manufacturer's Instructions For Use for specific weight and current limitations.

- For Endo FB3.0 and Endo SH2.0 devices, carefully insert and withdraw the device from trocar cannulas to avoid possible damage to the devices and/or injury to the patient.
- Use the device with an appropriate electrosurgical generator. See Table 1. Recommended Generator Settings.
- Surgery should be performed by persons with adequate training and preparation. Personnel should fully understand the nature and use of RF energy before performing electrosurgical procedures to avoid the risks of shock and burn hazards to both the patient and the operator and damage to the instrumentation.
- The cable on the device should be positioned in a way to avoid contact with the patient or other cables.
- Consult the operating and user manuals for other light sources, electrosurgical units, saline delivery systems and other ancillary devices for operating instructions, warnings and cautions prior to their use in the same surgical field as Medtronic's device.
- It is recommended that physicians utilize pre-clinical training, review of pertinent literature and other appropriate educational tools before attempting newer surgical procedures, such as endoscopic, laparoscopic, or thoracoscopic procedures.
- Examine the shipping carton, packaging, sterile barrier and device for any signs of transit damage. If there are any shortages, breakage or apparent damage, do not use the device. Return the device to Medtronic Advanced Energy LLC and use a new device.
- This device contains Phthalates.
- Read the instructions, warnings, and precautions provided with the electrosurgical device before using. Specific instructions for the generator are not included in this manual.
- Use the device with caution in the presence of pacemakers, as electrosurgical devices may cause interference with pacemakers or other active implants.
- The depth of effect is deeper and increases with time if the electrode is held stationary with less depth of effect if the electrode is moved over the tissue.
- High power settings result in deeper tissue effect than lower power settings.
- DO NOT use electrosurgery in the presence of flammable anesthetics or other flammable gases, near flammable fluids or objects, or in the presence of oxidizing agents, as fire could result.
- The patient should not come into contact with metal parts which are earthed or have an appreciable capacitance to earth (e.g. operating table supports, etc.).
- Skin to Skin contact (for example between the arms and body of the patient) should be avoided, for example by insertion of dry gauze.
- For surgical procedures where the RF current could flow through parts of the body having a relatively small cross sectional area, the use of bipolar techniques may be desirable in order to avoid unwanted tissue damage.
- Before using the device, confirm the following:
 - The cable from the device is connected to the monopolar output of the electrosurgical generator.
 - The dispersive electrode return pad is in place and appropriately connected, according to the manufacturer's instructions.
 - All electrical connections are tight, clean and dry.
 - All fluid connections are secure.
 - The generator is set at the desired power level.
 - The saline delivery tubing and device have been fully primed with sterile saline (0.9% NaCl) solution.

Setting Up the Monopolar Device

Warning:

- **Electric Shock Hazard.** Ensure that the device plug is correctly connected and that no metal pins are exposed.
- Use the lowest setting possible to achieve the desired tissue effect to avoid overtreatment, which could result in swelling, fluid, seroma or unintended tissue necrosis.
- Place any monitoring electrodes being used as far away as possible from the device.
 - Avoid needle monitoring electrodes.
 - Use monitoring systems incorporating high frequency current limiting devices.

Table 1. Recommended Generator Settings

Generator Make & Model	Mode/Setting
ConMed Excalibur, Sabre 2400	Standard Coag
ConMed System 5000	Standard Coag
Valleylab Force FX	Coag Fulgurate
Valleylab Force FX-C	Coag Fulgurate
Valleylab Force 2	Coag
Valleylab Force 30	Coag Fulgurate
Valleylab Force 40	Coag Fulgurate

- STEP 1:** Place dispersive electrode pad on patient, according to its manufacturer's instructions, and connect pad to appropriate generator.
- STEP 2:** Using aseptic technique, open the package and place contents on sterile field.
- STEP 3:** Set up and prime the saline delivery system according to manufacturer's Instructions For Use. Use an appropriate volume I.V. bag of sterile saline (0.9% NaCl) and a primary I.V. Tubing Set with a drip chamber.
- STEP 4:** Using aseptic technique, remove the vented luer cap from the device tubing. Connect tubing from the device to tubing from the saline delivery system, making sure the luer connections are secure.
- STEP 5:** Start saline flow and allow passage of saline through the tubing and device.
 - **Note:** In order to reduce set-up time, the device and tubing can be primed with a 60 ml syringe filled with sterile saline (0.9% NaCl). After priming the tubing, using aseptic technique, connect the saline delivery system line to the saline line provided with the device. Assume that the luer connection is secure.
- STEP 6:** Wait for saline to fill the device, and then invert the device so that the distal end faces the floor (tip end down).
- STEP 7:** Check flow of saline. Adjust the saline flow using the flow control on the I.V. set, or if using an infusion pump, set it to deliver 240 ml/hr. The saline should drip regularly from the end of device at approximately one to two drops per second.
- STEP 8:** Stop or pause the saline delivery until the device is ready to be used, using the pinch clamp provided with the device.
- STEP 9:** Connect the device cable to the monopolar receptacle on the generator. (Reference generator manufacturer's instructions.)
 - **Note:** For some generators, activation of the foot pedal will also activate coagulation. To use footswitch activation, some generators require that the device be connected to a footswitch controllable output. Some generators, particularly older ones, may not allow activation of the device with the footswitch at all. See Table 1. Recommended Generator Settings.

Set RF Power

Set generator to desired power setting. Start with the lowest recommended power setting and increase as needed to achieve the desired effect.

Activating and Operating the Monopolar Device

- STEP 1: Initiate saline flow and verify that saline is flowing from the tip of the device.
- STEP 2: Place tip of device on the tissue to be treated.
- STEP 3: Activate device by depressing button or foot pedal (see previous note) continuously for length of treatment, release button or foot pedal to cease treatment.
- STEP 4: Repeat as necessary by positioning the device over next area to be treated.
- STEP 5: Ensure that the tip is only in contact with tissue to be treated.
- STEP 6: For optimum performance, the tip must be kept free of debris.
- STEP 7: If saline flow stops during the electrosurgical procedure, stop using the device and attempt to resume saline flow. Ensure that the saline source is adequate and the saline delivery system is functioning properly. **If unable to resume saline flow, discontinue use and obtain a new device and return the used device to Medtronic Advanced Energy LLC.**

Disposing of the Monopolar Device

Dispose of the disposable device & the used saline bag according to procedures for your institution.

Precaution:

- The device and the saline bag will contain unused saline following use of the device. Take precautions to prevent the unused saline from flowing onto OR surfaces by placing handpiece into waste receptacle prior to opening pump head & removing the pump segment portion of the device.

Technical Description

General Description:

Monopolar electrosurgical device.
Employs RF energy and saline for hemostatic sealing, coagulation and, with the exception of the Endo FB3.0, blunt dissection of soft tissue.
See Table 1. Recommended Generator Settings.

General Information:

Sterile, EO
Single-use, Do Not Reuse
Non-Pyrogenic
Does not contain natural rubber latex
Caution: Read Instructions For Use (IFU) before using this device.

Physical Description:

Width: 1.05" (2.26 cm)
Length (without cables): 10" or 18.5" (25 cm or 47 cm)
Length of Electrical Cable: approximately 10 feet (3 meters)
Length of Saline Cable: approximately 16 feet (5 meters)
Weight (with cables): 4 oz. or 4.3 oz. (113 grams or 121 grams)

Operating Conditions:

Temperature: 50 °F to 86 °F (10 °C to 30 °C)
Humidity: 15% – 85%, non-condensing

Storage Conditions:

Temperature: 50 °F to 104 °F (10 °C to 40 °C)
Humidity: 15% – 85%, non-condensing

Handling Conditions:

Temperature: -20 °F to 122 °F (-29°C to 50 °C)
Humidity: 15% – 85%, non-condensing

Disposal:

Dispose of used devices as a biohazard.

Limited Express Warranty

IF, PRIOR TO THE DEVICE EXPIRATION DATE, THE STERILIZED DEVICE IS FOUND TO BE INOPERABLE DURING NORMAL AND PROPER USE IN ACCORDANCE WITH APPLICABLE INSTRUCTIONS, MEDTRONIC ADVANCED ENERGY LLC WILL REPLACE THE DEVICE AT NO CHARGE. MEDTRONIC ADVANCED ENERGY LLC MAKES NO OTHER WARRANTIES WITH RESPECT TO THE DEVICES AND EXPRESSLY DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, AS TO MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER MATTER. IN NO EVENT SHALL MEDTRONIC ADVANCED ENERGY LLC BE LIABLE FOR ANY CONSEQUENTIAL DAMAGES.

This Limited Express Warranty does not extend to a device where the user has compromised the sterile integrity of the package or the device is used after the device expiration date. A sterile package that has been compromised should not be used.

Customer Service

Please call Medtronic Customer Service at 866.777.9400 if you have any device returns or questions about Medtronic's device.

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